BVM EXT Tube User Manual

MANUFACTURED BY:
Micro BVM Systems Ltd.
11 Ruzin St. Jerusalem 93870 Israel

+972-2-6517191

info@microbvm.com

www.microbvm.com

(01) 07290016824186 1 (17) 250430 (10) 2001

EN

BVM EXT TUBE USER MANUAL

U.S. Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.

INTENDED USE

The BVM EXT Tube is intended for operation with standard manual pulmonary resuscitators (not included in this kit) to enable the use of resuscitator at a distance of up to 75 cm (30 inches) from the resuscitated patient.

The BVM EXT Tube is for single patient use and is not sterile. Its lifetime is 5 years.

The BVM EXT Tube is to be used by trained personnel only.

WARNINGS

- Oil or grease should not be used in close proximity to oxygen equipment; do not smoke or use open flames when oxygen is in use. Fire may result in both cases.
- Reduced or insufficient airflow may result in brain damage to the patient being ventilated.
- Improper function could lead to inadequate ventilation, resulting in serious injury or death.
- If efficient ventilation cannot be obtained with this device, the operator should immediately switch to mouth-to-mouth and/or mouth-to-nose resuscitation.
- Use in contaminated environments can be hazardous unless entrainment of contaminants is prevented.
- Oxygen flow should not exceed 35 liters per minute.

PRECAUTIONS AND CAUTIONS

- Only persons with adequate training should use the BVM EXT Tube.
- Users of this product must be familiar with the important information in the directors for use.
- Ensure to test each BVM resuscitator unit prior to use.
- This product is intended for single use only. Do not attempt to clean, decontaminate or sterilize.
- Never store the BVM EXT Tube in squeezed or compressed state, other than folded at the folding zones in the original packaging. Permanent deformation of the tube might occur, which may reduce the ventilation efficiency.

INSTRUCTIONS

1) Remove the BVM EXT Tube from its bag and expand the BVM EXT Tube from its compressed position to the required length by pulling the tube from its ends in opposite directions.

2) Attach the BVM EXT Tube to BVM resuscitator unit's front patient valve. Carry out all functional tests for the BVM unit while the BVM EXT Tube is attached.

3) Attach the patient's face mask or oropharyngeal airway to the BVM EXT Tube's patient valve housing.

4) Ventilate the patient. If the patient vomits immediately clear the BVM EXT Tube of vomitus. Use extreme caution if ventilating asthmatics or neonates.

FUNCTIONAL TESTS PRIOR TO USE

- Carry out functional tests for the BVM unit while the BVM EXT Tube is attached to it. Confirm that the BVM resuscitator unit passes all the tests.
- If there is an air leak from the BVM EXT Tube's patient valve, tighten the Expiratory connector and try again.
- If the BVM EXT Tube's patient valve is not fully functional, do not use it.

SPECIFICATIONS

Operating environmental limits	-18° C (0° F) to 50° C (122.9° F)
Storage environmental limits	-40° C (-40° F) to +60° C (140° F)
Gasses, other than air, that can be administered through the extension tube	Oxygen
Gas inlet tube connection	15 mm length x 6 mm O.D.
Expiratory resistance	3.4 cm H ₂ O is the resistance of the BVM unit
Delivered oxygen concentrations under various test conditions	55% (at 2 L/min) to 100% (at 8 L/min)
External dimensions of the extension tube (mm)	216.8 mm to 752.9 mm x 23.5 mm
Mass of the extension tube	61 grams

EXPLANATION OF SYMBOLS:

- Manufacturer
- European Authorized Representative
- Lot number
- Year of manufacture of active devices
- Limit temperature
- Medical device
- Use by

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